



Refrigerated Circulating Baths chiller

Available in 4 different models, the latest generation of Hally Refrigerated/Heated Circulating Baths uphold our reputation for innovation, product quality and refrigeration expertise. Our lab-and operator-friendly Circulating Baths set a new industry standard for design, performance and value.

Refrigerated Circulating Baths A versatile range of products to fit every need.

Capable of reaching temperatures as low as -35°C, Refrigerated/Heated Circulating Baths provide dependable temperature control for a wide variety of laboratory cooling applications.



Features :

- Working temperatures from -35° to +100°C
- Temperature stabilities to ±1°C
- Capacities from 7 to 50 liters
- Max flow rate 12 liter/ min.
- PID controller types
- Large bath openings and easy to clean surfaces
- Environmental Control cooling System



control

multi function PID level controller for temperature and soak time with delay time for compressor output are extremely easy to use. Clear navigational buttons and menu prompts guide you through setup, making operational and temperature set-point changes simple and precise.




pump

High speed submersible type pump compact technology reduces overall equipment noise, increases refrigeration efficiency and lowers energy consumption.



Energy Saving

Ultra cool technology carefully regulates the amount of cooling required, saving energy while delivering fast cool-downs and more precise control at elevated temperatures.



Safety

All models feature user-settable high-temperature safeties as well as redundant over-temperature protection.

Refrigerated/Heated Circulators come with the following accessories:

- Reservoir cover
- Bypass tubing
- Male inlet and outlet adapters for 3/8 tubing,

Responsible

An ISO 9001:2015 compliant Environmental Management System, recycling program and manufacturing are just some of the green initiatives Hally has undertaken to minimize our impact on the environment



CCL/QMS/23892/HI CE-2623

ISO 9001:2015 Certified Company

Note: Guarantee ceased in case of non use of suitable servo stabilizer & to be stored in A.C. room.
Due to continues R&D model/specification may be change.

HALLY INSTRUMENTS MUMBAI

Unit No. 112, 1st Flr., 'C' wing, Vasupooja Industrial Estate , Laxmi Nagar., Goregoan (w) Mumbai -400062. Tel. : 2876 9190, Fax : 2877 7670,
E-mail : hally_instruments@yahoo.co.in. Web : www.hally.in